



4141 MacArthur Blvd. • Newport Beach, CA 92660
 800-854-7256 • Fax 800-411-9722 • glidewell.com

Dr. Name _____ Phone # _____

Acct. # _____ Patient ID/Name _____
First Last

Address/Email _____

Deliver by 5 p.m. on _____

Enclosed with Case: Impressions Models Bite Photos Other: _____

SPECIAL OFFER*
 Save 25% on your next crown, bridge or veneer case (up to 3 units) by enclosing this Rx along with your physical impressions.
Offer expires May 31, 2024. *Price does not include shipping or applicable taxes. Limit one special offer per case. No copies accepted. Special pricing not valid with any other offer.

ZIRCONIA/ALL-CERAMIC RESTORATIONS

- BruxZir Full-Strength† (> 1,000 MPa)
- IPS e.max
- BruxZir Esthetic (870 MPa)
- Lava Crowns & Bridges
- Obsidian All-Ceramic
- Bilayered Clinical Zirconia

VENEERS

- Obsidian veneer†
- BruxZir Esthetic veneer
- IPS e.max veneer
- Layered IPS e.max veneer

PFM

- Obsidian to Non-Precious†
- Obsidian to White Noble
- Porcelain to OcclusalGold YHN
- Obsidian to White High Noble

FULL-CAST RESTORATIONS

- Noble-Cast 45 YN (40% Au)
- White Noble
- Noble-Cast 60 YHN (57.5% Au)†
- White High Noble (40% Au)
- Non-Precious

FINAL SHADE



STUMP SHADE

Must indicate prepped tooth shade for all-ceramics

OCCUSAL STAINING

- None
- Light†
- Medium
- Dark

PONTIC DESIGN



MARGIN AND METAL DESIGN




IF NO OCCLUSAL CLEARANCE

- Call doctor
- Spot opposing
- Metal occlusion
- Metal island
- Make this a permanent note

WEB Rx

25% off crown & bridge restoration



Signature _____

License # _____ Date _____

Submission of this Rx constitutes agreement with limited warranty terms and conditions.
 For details, visit glidewell.com/policies-and-warranties.

TERMS AND WARRANTY INFORMATION



All Restorations Made in the USA

We honor VISA, MASTERCARD, AMEX and DISCOVER.




TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: Glidewell is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit glidewell.com/policies-and-warranties.

PREPARATION GUIDELINES



BruxZir Esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

BruxZir Full-Strength

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

THE JIM GLIDEWELL SIGNATURE ANTERIOR SERIES

Choose the desired esthetic outcome for your patient

SQUARE-TAPERED



Exemplify boldness and youthfulness.

SOFT-SQUARED



Show energetic professionalism.

OVOID



Convey charm and softness.

TRAPEZOID



Demonstrate confidence.

SQUARED



Strike an athletic tone.

TRIANGLE-TAPERED



Exhibit mature experience.