



Glidewell HT™
IMPLANT SYSTEM

Product Catalog

glidewell-ht.com

5	Technical Specifications	22	Drivers & Accessories
6	Implants	23	Guided Surgical Instrumentation Kit
7	Healing Components	24	Guided Surgical Kit Overview
8	Impression Components	26	Guided Drills
10	Stock Abutments	28	Guided Taps & Accessories
13	Multi-Unit Abutments	30	Prosthetic Instrumentation Kit
15	Multi-Unit Accessories	31	Prosthetic Instrumentation Kit Overview
16	LOCATOR® Abutment System	32	Product Packaging
18	Standard Surgical Instrumentation Kits	33	Policies and Warranty
19	Surgical Kit Overview	34	Parts Index
20	Drills & Taps		

Glidewell HT™ Implant System is a trademark of Prismatic Dentalcraft, Inc.

Copyright © 2025, Prismatic Dentalcraft, Inc. Prismatic Dentalcraft, Inc. is not responsible for any damages or other liabilities (including attorney fees) resulting, or claimed to result in whole or in part, from actual or alleged problems arising out of the use of this information. The techniques, procedures and theories presented herein are provided in good faith and believed to be correct as of the date hereof. Any dental professional viewing this presentation must make his or her own decisions about the use of the materials and techniques for specific situations.

No representations as to the completeness or accuracy of this information is given, and no representations or warranties, either expressed or implied, of merchantability, fitness for a particular purpose or of any other nature are made here under with respect to the information or the product to which information refers.

All third-party trademarks are property of their respective owners.

The Glidewell HT™ Implant System

The **Glidewell HT™ Implant System**, formerly known as the Hahn™ Tapered Implant System, is a premium implant solution that simplifies surgery and provides you with unrivaled support from the most experienced dental lab in the U.S. With a 99.2% long-term success rate backed by an implant-to-crown lifetime warranty, the Glidewell HT Implant promotes success in implant dentistry while lowering your surgical and restorative costs.¹

- **Simple and efficient** — Easy to use, with a streamlined surgical protocol, length-specific drills, and new-and-improved surgical kits based on eight years of clinical feedback
- **Cut your costs** — Priced at a fraction of comparable implant systems — and saves you 20% on your lab bill when you restore your implant case with Glidewell*
- **Clinically proven** — 99.2% success rate and 0.2 mm mean bone loss¹
- **High primary stability** — Deep, sharp threads maximize initial stability and engage bone where directed

“The **Glidewell HT™ Implant System** is engineered to help dentists provide implant treatment for more patients through **ease of use, reduced costs** and our **unwavering commitment to support your practice** — from implant placement to final restoration.”

Jim Glidewell, CDT
Founder and President of Glidewell



About the Manufacturer

Prismatik Dentalcraft was established in 2006 and includes a carefully assembled team of experts with a proven track record in the design, engineering and manufacture of dental implants. Bolstered by a support staff of highly respected researchers, material scientists, clinical specialists, and dental technicians, Prismatik is dedicated to advancing implant therapies by combining proven treatment protocols with progressive materials, technologies and techniques.



Vertical Integration

Our ownership of the entire manufacturing process behind our implant products ensures quality and helps reduce costs for our customers.



State-of-the-Art Equipment

Our Swiss-type lathes and CNC milling machines are ideal for implants and prosthetics requiring extreme precision.



Made in the USA

Our ISO-certified facility in Irvine, California, operates under FDA Current Good Manufacturing Practices (CGMPs).

*Discount offered only at Glidewell and cannot be combined with any other special offers. Case must include an implant-level or multi-unit abutment-level impression with a Glidewell HT transfer coping or a digital scan with a Glidewell HT scan body. Impressions over cementable abutments are not eligible for discount.

1. Kerr M, Allen B, Park N. Clinical and radiographic evaluation of tapered implants with an aggressive reverse buttress thread and crestal microthreads: a retrospective study. For the full report, visit glidewell.com/ht-2-year.

PROVEN, TIME-TESTED IMPLANT DESIGN FEATURES

99.2%

SUCCESS RATE¹

0.2 mm

MEAN BONE LOSS IN 2-YEAR STUDY¹

REDUCED MARGINAL BONE STRESS²

Convergent implant neck design minimizes forces on cortical bone.

ENHANCED BIOLOGICAL SEAL^{2,13}

Platform-shifted abutment interface and conical connection provide excellent seal, stability and tissue preservation.

RBM SURFACE TECHNOLOGY³⁻⁵

Resorbable blast media (RBM) surface treatment results in textured, ultra-clean surface that promotes osseointegration.

SWIFT INSERTION⁶

Dual-lead thread pattern with self-tapping grooves.



CRESTAL BONE AND SOFT-TISSUE PRESERVATION^{2,7-9}

Machined collar surface leads to 25-times less biofilm attachment and mitigates peri-implantitis. Coronal microthreads reduce crestal bone loss.



HIGH PRIMARY STABILITY AND DIRECTIONAL CONTROL DURING PLACEMENT^{10,11}

Deep, sharp threads maximize initial stability for immediate load opportunities and long-term success.


TAPERED IMPLANT BODY¹²

Tapered body shape — a design feature originally conceived by Dr. Jack Hahn — eases placement in anatomically constricted areas.

Technical Specifications



Glidewell HT™ Implants allow for precise control during placement, engage a maximum amount of bone, and achieve a high degree of primary stability in a wide variety of clinical situations. The unique thread pattern with self-tapping grooves is sharply angled for swift insertion. All Glidewell HT Implants feature a conical internal hex connection for a secure prosthetic seal. Matching components are color-coded by implant diameter and/or restorative platform. To help minimize the overall number of system components, 3.5 mm and 4.3 mm diameter implants share the same restorative platform.




Description	Code
Glidewell HT™ Implant	
Ø3.0 x 11.5 mm	70-1189-IMP0001
Ø3.0 x 13 mm	70-1189-IMP0002
Ø3.0 x 16 mm	70-1189-IMP0003
Ø3.5 x 8 mm	70-1189-IMP0004
Ø3.5 x 10 mm	70-1189-IMP0005
Ø3.5 x 11.5 mm	70-1189-IMP0006
Ø3.5 x 13 mm	70-1189-IMP0007
Ø3.5 x 16 mm	70-1189-IMP0008
Ø4.3 x 8 mm	70-1189-IMP0009
Ø4.3 x 10 mm	70-1189-IMP0010
Ø4.3 x 11.5 mm	70-1189-IMP0011
Ø4.3 x 13 mm	70-1189-IMP0012
Ø4.3 x 16 mm	70-1189-IMP0013
Ø5.0 x 8 mm	70-1189-IMP0014
Ø5.0 x 10 mm	70-1189-IMP0015
Ø5.0 x 11.5 mm	70-1189-IMP0016
Ø5.0 x 13 mm	70-1189-IMP0017
Ø5.0 x 16 mm	70-1189-IMP0018
Ø6.0 x 8 mm	70-1189-IMP0022
Ø6.0 x 10 mm	70-1189-IMP0023
Ø6.0 x 11.5 mm	70-1189-IMP0024
Ø7.0 x 8 mm	70-1189-IMP0019
Ø7.0 x 10 mm	70-1189-IMP0020
Ø7.0 x 11.5 mm	70-1189-IMP0021

(Implants include Titanium Cover Screw.)

Description

Code

Glidewell HT™ Implant Healing Abutment



3 mmH - Ø3.0 Implant	70-1190-IMP0001
5 mmH - Ø3.0 Implant	70-1190-IMP0002
7 mmH - Ø3.0 Implant	70-1190-IMP0003
3 mmH - Ø3.5 Implant	70-1190-IMP0004
5 mmH - Ø3.5 Implant	70-1190-IMP0005
7 mmH - Ø3.5 Implant	70-1190-IMP0006
3 mmH - Ø4.3 Implant	70-1190-IMP0007
5 mmH - Ø4.3 Implant	70-1190-IMP0008
7 mmH - Ø4.3 Implant	70-1190-IMP0009
3 mmH - Ø5.0 Implant	70-1190-IMP0010
5 mmH - Ø5.0 Implant	70-1190-IMP0011
7 mmH - Ø5.0 Implant	70-1190-IMP0012
3 mmH - Ø6.0 Implant	70-1190-IMP0015
5 mmH - Ø6.0 Implant	70-1190-IMP0016
3 mmH - Ø7.0 Implant	70-1190-IMP0013
5 mmH - Ø7.0 Implant	70-1190-IMP0014

(Healing Abutments require a Glidewell HT™ Implant Prosthetic Driver for placement; please see page 22.)

Glidewell HT™ Implant Temporary Abutment



Ø3.0 Implant	70-1190-PRA0025
Ø3.5/4.3 Implant	70-1190-PRA0026
Ø5.0 Implant	70-1190-PRA0027
Ø6.0 Implant	70-1190-PRA0062
Ø7.0 Implant	70-1190-PRA0028

(Temporary Abutments include Titanium Screw.)

Description

Code


Glidewell HT™ Implant Non-Engaging Temporary Abutment



Ø3.0 Implant	70-1190-PRA0021
Ø3.5/4.3 Implant	70-1190-PRA0022
Ø5.0 Implant	70-1190-PRA0023
Ø6.0 Implant	70-1190-PRA0063
Ø7.0 Implant	70-1190-PRA0024

(Temporary Abutments include Titanium Screw.)

Glidewell HT™ Implant Closed-Tray Impression Coping



3 mmH - Ø3.0 Implant	70-1190-PRC0001
5 mmH - Ø3.0 Implant	70-1190-PRC0002
7 mmH - Ø3.0 Implant	70-1190-PRC0003
3 mmH - Ø3.5 Implant	70-1190-PRC0004
5 mmH - Ø3.5 Implant	70-1190-PRC0005
7 mmH - Ø3.5 Implant	70-1190-PRC0006
3 mmH - Ø4.3 Implant	70-1190-PRC0007
5 mmH - Ø4.3 Implant	70-1190-PRC0008
7 mmH - Ø4.3 Implant	70-1190-PRC0009
3 mmH - Ø5.0 Implant	70-1190-PRC0010
5 mmH - Ø5.0 Implant	70-1190-PRC0011
7 mmH - Ø5.0 Implant	70-1190-PRC0012
3 mmH - Ø6.0 Implant	70-1190-PRC0058
5 mmH - Ø6.0 Implant	70-1190-PRC0059
3 mmH - Ø7.0 Implant	70-1190-PRC0013
5 mmH - Ø7.0 Implant	70-1190-PRC0014

(Closed-Tray Impression Copings include Closed-Tray Screw.)


Description

Code


Glidewell HT™ Implant Open-Tray Impression Coping




3 mmH - Ø3.0 Implant	70-1190-PRC0015
5 mmH - Ø3.0 Implant	70-1190-PRC0016
7 mmH - Ø3.0 Implant	70-1190-PRC0017




3 mmH - Ø3.5 Implant	70-1190-PRC0018
5 mmH - Ø3.5 Implant	70-1190-PRC0019
7 mmH - Ø3.5 Implant	70-1190-PRC0020



3 mmH - Ø4.3 Implant	70-1190-PRC0021
5 mmH - Ø4.3 Implant	70-1190-PRC0022
7 mmH - Ø4.3 Implant	70-1190-PRC0023



3 mmH - Ø5.0 Implant	70-1190-PRC0024
5 mmH - Ø5.0 Implant	70-1190-PRC0025
7 mmH - Ø5.0 Implant	70-1190-PRC0026



3 mmH - Ø6.0 Implant	70-1190-PRC0060
5 mmH - Ø6.0 Implant	70-1190-PRC0061



3 mmH - Ø7.0 Implant	70-1190-PRC0027
5 mmH - Ø7.0 Implant	70-1190-PRC0028

(Open-Tray Impression Copings include Open-Tray Screw.)

Impression Components

Description		Code
Glidewell HT™ Implant Titanium Scan Body		
	Ø3.0 Implant	70-1190-PRA0054
	Ø3.5/4.3 Implant	70-1190-PRA0055
	Ø5.0 Implant	70-1190-PRA0056
	Ø6.0 Implant	70-1190-PRA0064
	Ø7.0 Implant	70-1190-PRA0057

Glidewell HT™ Implant Analog

	Ø3.0 Implant	70-1190-PRC0035
	Ø3.5 Implant	70-1190-PRC0036
	Ø4.3 Implant	70-1190-PRC0037
	Ø5.0 Implant	70-1190-PRC0038
	Ø6.0 Implant	70-1190-PRC0062
	Ø7.0 Implant	70-1190-PRC0039

Stock Abutments

Glidewell HT™ Implant Anterior Esthetic Abutment

	Ø3.0 Implant	70-1190-PRA0037
	Ø3.5/4.3 Implant	70-1190-PRA0038
	Ø5.0 Implant	70-1190-PRA0039

Glidewell HT™ Implant Posterior Esthetic Abutment

	Ø5.0 Implant	70-1190-PRA0040
	Ø6.0 Implant	70-1190-PRA0075
	Ø7.0 Implant	70-1190-PRA0041

(Esthetic Abutments include Titanium Screw.)

Description

Code

Glidewell HT™ Implant 15° Anterior Esthetic Abutment



Ø3.0 Implant	70-1190-PRA0042
Ø3.5/4.3 Implant	70-1190-PRA0043
Ø5.0 Implant	70-1190-PRA0044

Glidewell HT™ Implant 15° Posterior Esthetic Abutment



Ø5.0 Implant	70-1190-PRA0045
--------------	-----------------

(Esthetic Abutments include Titanium Screw.)

Glidewell HT™ Implant Non-Engaging Titanium Abutment



4.5 mmH - Ø3.0 Implant	70-1190-PRA0001
4.5 mmH - Ø3.5/4.3 Implant	70-1190-PRA0002
4.5 mmH - Ø5.0 Implant	70-1190-PRA0003



4.5 mmH - Ø6.0 Implant	70-1190-PRA0065
4.5 mmH - Ø7.0 Implant	70-1190-PRA0004

Glidewell HT™ Implant Titanium Abutment



4.5 mmH - Ø3.0 Implant	70-1190-PRA0005
4.5 mmH - Ø3.5/4.3 Implant	70-1190-PRA0006
4.5 mmH - Ø5.0 Implant	70-1190-PRA0007



4.5 mmH - Ø6.0 Implant	70-1190-PRA0066
4.5 mmH - Ø7.0 Implant	70-1190-PRA0008



6 mmH - Ø3.0 Implant	70-1190-PRA0009
6 mmH - Ø3.5/4.3 Implant	70-1190-PRA0010
6 mmH - Ø5.0 Implant	70-1190-PRA0011



6 mmH - Ø6.0 Implant	70-1190-PRA0067
6 mmH - Ø7.0 Implant	70-1190-PRA0012

(Titanium Screw sold separately; please see page 13.)

Description		Code
Glidewell HT™ Implant UCLA Gold		
	Ø3.0 Implant	70-1190-PRA0046
	Ø3.5/4.3 Implant	70-1190-PRA0047
	Ø5.0 Implant	70-1190-PRA0048
	Ø6.0 Implant	70-1190-PRA0070
	Ø7.0 Implant	70-1190-PRA0049
Glidewell HT™ Implant Non-Engaging UCLA Gold		
	Ø3.0 Implant	70-1190-PRA0050
	Ø3.5/4.3 Implant	70-1190-PRA0051
	Ø5.0 Implant	70-1190-PRA0052
	Ø6.0 Implant	70-1190-PRA0073
	Ø7.0 Implant	70-1190-PRA0053
Glidewell HT™ Implant UCLA Plastic		
	Ø3.0 Implant	70-1190-PRA0029
	Ø3.5/4.3 Implant	70-1190-PRA0030
	Ø5.0 Implant	70-1190-PRA0031
	Ø6.0 Implant	70-1190-PRA0074
	Ø7.0 Implant	70-1190-PRA0032
Glidewell HT™ Implant Non-Engaging UCLA Plastic		
	Ø3.0 Implant	70-1190-PRA0033
	Ø3.5/4.3 Implant	70-1190-PRA0034
	Ø5.0 Implant	70-1190-PRA0035
	Ø6.0 Implant	70-1190-PRA0071
	Ø7.0 Implant	70-1190-PRA0036

(UCLA Abutments include Titanium Screw.)

Stock Abutments

	Description	Code
Glidewell HT™ Implant Titanium Screw and Guide Pin		
 70-1190-PRC0029	Ø3.0 Implant	70-1190-PRC0029
 70-1190-PRC0030	Ø3.5/4.3 Implant	70-1190-PRC0030
 70-1190-PRC0064 70-1190-PRC0031	Ø5.0/6.0 Implant Ø5.0/7.0 Implant	70-1190-PRC0064 70-1190-PRC0031
 70-1190-PRC0047	Ø3.0 Implant	70-1190-PRC0047
 70-1190-PRC0048	Ø3.5/4.3 Implant	70-1190-PRC0048
 70-1190-PRC0065 70-1190-PRC0049	Ø5.0/6.0 Implant Ø5.0/7.0 Implant	70-1190-PRC0065 70-1190-PRC0049

Multi-Unit Abutments

Glidewell HT™ Implant Multi-Unit Abutment		
 70-1190-PRS0011	2.5 mmH - Ø3.0 Implant	70-1190-PRS0011
 70-1190-PRS0014	3.5 mmH - Ø3.0 Implant	70-1190-PRS0014
 70-1190-PRS0009	1.5 mmH - Ø3.5/4.3 Implant	70-1190-PRS0009
 70-1190-PRS0012	2.5 mmH - Ø3.5/4.3 Implant	70-1190-PRS0012
 70-1190-PRS0015	3.5 mmH - Ø3.5/4.3 Implant	70-1190-PRS0015
 70-1190-PRS0010	1.5 mmH - Ø5.0 Implant	70-1190-PRS0010
 70-1190-PRS0013	2.5 mmH - Ø5.0 Implant	70-1190-PRS0013
 70-1190-PRS0016	3.5 mmH - Ø5.0 Implant	70-1190-PRS0016
 70-1190-PRS0017	4.5 mmH - Ø5.0 Implant	70-1190-PRS0017
 70-1190-PRS0019	1.5 mmH - Ø6.0 Implant	70-1190-PRS0019
 70-1190-PRS0020	2.5 mmH - Ø6.0 Implant	70-1190-PRS0020
 70-1190-PRS0021	3.5 mmH - Ø6.0 Implant	70-1190-PRS0021
 70-1190-PRS0022	4.5 mmH - Ø6.0 Implant	70-1190-PRS0022

Description		Code
Glidewell HT™ Implant 17° Multi-Unit Abutment		
 70-1190-PRS0001	2.5 mmH - Ø3.5/4.3 Implant	70-1190-PRS0001
 70-1190-PRS0003	3.5 mmH - Ø3.5/4.3 Implant	70-1190-PRS0003
 70-1190-PRS0002	2.5 mmH - Ø5.0 Implant	70-1190-PRS0002
 70-1190-PRS0004	3.5 mmH - Ø5.0 Implant	70-1190-PRS0004

Glidewell HT™ Implant 30° Multi-Unit Abutment		
 70-1190-PRS0005	3.5 mmH - Ø3.5/4.3 Implant	70-1190-PRS0005
 70-1190-PRS0007	4.5 mmH - Ø3.5/4.3 Implant	70-1190-PRS0007
 70-1190-PRS0006	3.5 mmH - Ø5.0 Implant	70-1190-PRS0006
 70-1190-PRS0008	4.5 mmH - Ø5.0 Implant	70-1190-PRS0008

(Angled Multi-Unit Abutments include Angled Abutment Screw. Prosthetic Screw sold separately; please see page 15.)



Glidewell HT™ Implant Multi-Unit Coping

Multi-Unit Coping 4.5mmH, Non-Engaging	70-1190-PRS0018
--	-----------------



Glidewell HT™ Implant Angled Multi-Unit Abutment Screw

Ø3.5/4.3 Implant	70-1190-PRC0032
Ø5.0 Implant	70-1190-PRC0033






Multi-Unit Accessories

The Glidewell HT™ Implant System proudly features Inclusive® and Newport Surgical™ brand multi-unit accessories.

	Description	Code
Glidewell HT™ Implant Multi-Unit		
 70-1190-PRC0046	Titanium Scan Body with Prosthetic Screw	70-1190-PRC0046
Inclusive® Multi-Unit		
 70-1074-PRS0044	Temporary Healing Cap with Prosthetic Screw - Unigrip™ Driver	70-1074-PRS0044
 70-1074-PRS0057	Titanium Temporary with Prosthetic Screw - Unigrip™ Driver	70-1074-PRS0057
 70-1072-PRC0043	Prosthetic Screw - Unigrip™ Driver	70-1072-PRC0043
 70-1074-PRS0030	Gold/Plastic Coping	70-1074-PRS0030
 70-1074-PRS0018	Gold Bar Coping	70-1074-PRS0018
 70-1074-PRS0069	Plastic Coping with Prosthetic Screw - Unigrip™ Driver	70-1074-PRS0069
 70-1074-PRS0022	Impression Coping, Open-Tray, with Guide Pin - Unigrip™ Driver	70-1074-PRS0022
 70-1074-PRS0021	Impression Coping, Closed-Tray	70-1074-PRS0021
 70-1074-PRS0016	Abutment Analog	70-1074-PRS0016
 70-1074-PRS0019	Guide Pin - Unigrip™ Driver	70-1074-PRS0019
<i>Unigrip is a trademark of Nobel Biocare Services AG.</i>		
Newport Surgical™ Multi-Unit		
 70-1071-SRG0234	Newport Surgical Multi-Unit Driver, 4mm Square	70-1071-SRG0234


LOCATOR® Abutment System

The Glidewell HT™ Implant System proudly features Locator® brand overdenture attachments.

Description	Code
LOCATOR® Abutment	
 70-1075-ZST0065	Conical Connection Ø3.0 x 1.0 mmH Cuff [#2212] 70-1075-ZST0065
 70-1075-ZST0074	Conical Connection Ø3.0 x 2.0 mmH Cuff [#2213] 70-1075-ZST0074
 70-1075-ZST0066	Conical Connection Ø3.0 x 3.0 mmH Cuff [#2214] 70-1075-ZST0066
 70-1075-ZST0032	Conical Connection Ø3.5/4.3 x 1.0 mmH Cuff [#2069] 70-1075-ZST0032
 70-1075-ZST0034	Conical Connection Ø3.5/4.3 x 3.0 mmH Cuff [#2071] 70-1075-ZST0034
 70-1075-ZST0036	Conical Connection Ø3.5/4.3 x 5.0 mmH Cuff [#2073] 70-1075-ZST0036
 70-1075-ZST0071	Conical Connection Ø5.0 x 1.0 mmH Cuff [#2075] 70-1075-ZST0071
 70-1075-ZST0038	Conical Connection Ø5.0 x 3.0 mmH Cuff [#2077] 70-1075-ZST0038
 70-1075-ZST0040	Conical Connection Ø5.0 x 5.0 mmH Cuff [#2079] 70-1075-ZST0040

LOCATOR Abutment System

The Glidewell HT Implant System proudly features Locator brand overdenture attachments.

	Description	Code
LOCATOR®		
 70-1075-ZST0024	Male Processing Package (2-pack) [#8519-2] <i>(Includes: Denture Cap with Black Processing Male, Clear Nylon Male (5 lbs), Pink (LR) Less-Retentive Nylon Male (3 lbs), Blue Extra-Light-Retention Nylon Male (1.5 lbs), and White Block-out Spacer.)</i>	70-1075-ZST0024
 70-1075-ZST0025	Extended Range Male Processing Package (2-pack) [#8540-2] <i>(Includes: Denture Cap with Black Processing Male, Extended-Range Replacement Males: Green (3-4 lbs), Orange (2 lbs), Red (0.5-1.5 lbs), and White Block-out Spacer.)</i>	70-1075-ZST0025
 70-1075-ZST0007	Impression Coping (4-pack) [#8505]	70-1075-ZST0007
 70-1075-ZST0009	Female Analog 4 mmD (4-pack) [#8530]	70-1075-ZST0009
 70-1075-ZST0005	Core Tool Abutment Holder Sleeve (4-pack) [#8394]	70-1075-ZST0005
 70-1075-ZST0026	Angle Measurement Guide [#9530]	70-1075-ZST0026
 70-1075-ZST0010	Processing Spacer (4-pack) [#8569]	70-1075-ZST0010
 70-1075-ZST0004	Core Tool (Complete) [#8393]	70-1075-ZST0004
 70-1075-ZST0006	Core Tool Male Removal (Tip End) [#8397]	70-1075-ZST0006
 70-1075-ZST0008	Parallel Post (4-pack) [#8517]	70-1075-ZST0008
 70-1075-ZST0003	Core Tool Abutment Driver (Gold End) [#8390]	70-1075-ZST0003
 70-1075-ZST0001	35 Ncm Torque Wrench Insert Driver (15 mm) [#8260]	70-1075-ZST0001
 70-1075-ZST0002	35 Ncm Torque Wrench Insert Driver (21 mm) [#8280]	70-1075-ZST0002
 70-1075-ZST0020	Latch Type Torque Wrench Insert Driver (23 mm) [#8913]	70-1075-ZST0020
 70-1075-ZST0021	Latch Type Torque Wrench Insert Driver (29 mm) [#8914]	70-1075-ZST0021
 70-1075-ZST0022	Square Drive Torque Wrench Insert Driver (15 mm) [#8926]	70-1075-ZST0022
 70-1075-ZST0023	Square Drive Torque Wrench Insert Driver (21 mm) [#8927]	70-1075-ZST0023

LOCATOR is a registered trademark of Zest Anchors, LLC.

Standard Surgical Instrumentation Kits

Features and Benefits

The instrumentation kits allow the clinician to easily organize, store and transport the tooling of the Glidewell HT™ Implant System. The autoclavable kits have been designed for ease of use, with component labels for identification.

The Glidewell HT Implant Surgical Kit includes instrumentation that is machined from corrosion-resistant, surgical stainless steel, and features standard connectivity. The color-coded surgical kit clearly maps the drilling protocol for each implant. Drills are arranged from left to right in order of increasing diameter, with color-coded bands on the shaping drills to indicate the corresponding diameter of the Glidewell HT Implant. For added safety, each shaping drill is implant-specific to precisely determine both diameter and depth of the osteotomy. **Please note that Ø6.0 mm and Ø7.0 mm Shaping Drills, Screw Tap, Direction Indicator, and Driver are sold separately (refer to diagram, page 19).**


All instrumentation is manufactured in the U.S.A., Switzerland, Denmark or Germany. For specific country of origin, please refer to the individual product label.


Sterility

Please refer to the Glidewell HT Implant System Surgical Manual for cleaning and sterilization instructions.

Description	Code
Glidewell HT™ Implant Surgical Kit	70-1071-SRG0396







Glidewell HT™ Implant Twist Drill

Description	Code
Ø1.5 x 8 mm	70-1071-SRG0266
Ø2.4/1.5 mm	70-1071-SRG0264
Ø2.8/2.4 x 8 mm	70-1071-SRG0268
Ø2.8/2.4 x 10 mm	70-1071-SRG0269
Ø2.8/2.4 x 11.5 mm	70-1071-SRG0270
Ø2.8/2.4 x 13 mm	70-1071-SRG0271
Ø2.8/2.4 x 16 mm	70-1071-SRG0272



Glidewell HT™ Implant Shaping Drill

Ø3.5 x 8 mm	70-1071-SRG0273
Ø3.5 x 10 mm	70-1071-SRG0274
Ø3.5 x 11.5 mm	70-1071-SRG0275
Ø3.5 x 13 mm	70-1071-SRG0276
Ø3.5 x 16 mm	70-1071-SRG0277
Ø4.3 x 8 mm	70-1071-SRG0278
Ø4.3 x 10 mm	70-1071-SRG0279
Ø4.3 x 11.5 mm	70-1071-SRG0280
Ø4.3 x 13 mm	70-1071-SRG0281
Ø4.3 x 16 mm	70-1071-SRG0282
Ø5.0 x 8 mm	70-1071-SRG0283
Ø5.0 x 10 mm	70-1071-SRG0284



70-1071-SRG0285



70-1071-SRG0286



70-1071-SRG0287



70-1071-SRG0415



70-1071-SRG0416



70-1071-SRG0417



70-1071-SRG0288



70-1071-SRG0289



70-1071-SRG0290

Glidewell HT™ Implant Shaping Drill

Description	Code
Ø5.0 x 11.5 mm	70-1071-SRG0285
Ø5.0 x 13 mm	70-1071-SRG0286
Ø5.0 x 16 mm	70-1071-SRG0287
Ø6.0 x 8 mm	70-1071-SRG0415
Ø6.0 x 10 mm	70-1071-SRG0416
Ø6.0 x 11.5 mm	70-1071-SRG0417
Ø7.0 x 8 mm	70-1071-SRG0288
Ø7.0 x 10 mm	70-1071-SRG0289
Ø7.0 x 11.5 mm	70-1071-SRG0290



70-1071-SRG0291



70-1071-SRG0292



70-1071-SRG0293



70-1071-SRG0294




70-1071-SRG0414



70-1071-SRG0295


Glidewell HT™ Implant Screw Tap

Description	Code
Ø3.0 mm	70-1071-SRG0291
Ø3.5 mm	70-1071-SRG0292
Ø4.3 mm	70-1071-SRG0293
Ø5.0 mm	70-1071-SRG0294
Ø6.0 mm	70-1071-SRG0414
Ø7.0 mm	70-1071-SRG0295




Glidewell HT™ Implant Driver

Description	Code
Ø3.0 Long	70-1071-SRG0296
Ø3.5/4.3 Short	70-1071-SRG0297
Ø3.5/4.3 Long	70-1071-SRG0298
Ø5.0 Short	70-1071-SRG0299
Ø5.0 Long	70-1071-SRG0300
Ø6.0 Short	70-1071-SRG0419
Ø7.0 Short	70-1071-SRG0301




Glidewell HT™ Implant Direction Indicator

Ø3.0 mm	70-1071-SRG0308
Ø3.5 mm	70-1071-SRG0309
Ø4.3 mm	70-1071-SRG0310
Ø5.0 mm	70-1071-SRG0311
Ø6.0 mm	70-1071-SRG0413
Ø7.0 mm	70-1071-SRG0312



Glidewell HT™ Implant Surgical Adaptor - Short	70-1071-SRG0307
--	-----------------

Glidewell HT™ Torque Wrench with Driver Handle, Hex and Square	70-C35120
--	-----------



Inclusive® Drill Extender	70-1071-SRG0014
---------------------------	-----------------

Glidewell HT™ Prosthetic Driver - 30 mm	70-1190-PRC0056
---	-----------------

Glidewell HT™ Prosthetic Driver - 18 mm	70-1190-PRC0054
---	-----------------



Glidewell HT™ Implant Radiographic Template	70-1071-RAD0001
---	-----------------

Guided Surgical Instrumentation Kit

Features and Benefits

The Glidewell HT™ Implant Guided Surgery System allows clinicians to deliver Glidewell HT Implants with the highest level of predictability, accuracy, and efficiency. Designed for ease of use, including enhancements based on eight years of clinical feedback, this system works in conjunction with digital treatment planning and case-specific surgical guides.

The Glidewell HT™ Implant Guided Surgical Kit includes instrumentation that is machined from high-quality, corrosion-resistant, surgical stainless steel. Color-coded bands correspond to implant diameter, matching the color of the implant carrier. Tissue punches match the diameter of the prescribed implant. Drills are diameter-specific and feature a flange stop for depth control. For added safety, each shaping drill is implant-specific to precisely determine both diameter and depth of the osteotomy. Drills are arranged from left to right in order of increasing diameter, following the recommended drilling sequence.

All instrumentation is manufactured in the U.S.A., Switzerland, Denmark or Israel. For specific country of origin, please refer to the individual product label.



Sterility

Please refer to the Glidewell HT Implant System Surgical Manual for cleaning and sterilization instructions.


Description	Code
Glidewell HT™ Implant Guided Surgical Kit	70-1071-SRG0395



Guided Surgical Kit Overview



Guided Surgical Kit Overview





70-1071-SRG0313



70-1071-SRG0314



70-1071-SRG0315



70-1071-SRG0316

Description

Code

Glidewell HT™ Implant Guided Tissue Punch

Ø3.0 mm

70-1071-SRG0313

Ø3.5 mm

70-1071-SRG0314

Ø4.3 mm

70-1071-SRG0315

Ø5.0 mm

70-1071-SRG0316



70-1071-SRG0317



70-1071-SRG0318



70-1071-SRG0319

Glidewell HT™ Implant Guided Alignment Drill

Ø3.0 mm

70-1071-SRG0317

Ø3.5 mm

70-1071-SRG0318

Ø4.3/5.0 mm

70-1071-SRG0319



70-1071-SRG0387



70-1071-SRG0388



70-1071-SRG0389



70-1071-SRG0390



70-1071-SRG0391



70-1071-SRG0392



70-1071-SRG0393

Glidewell HT™ Implant Guided Pilot Drill

Ø3.0 x 8 mm

70-1071-SRG0387

Ø3.5 x 8 mm

70-1071-SRG0388

Ø3.5 x 10 mm

70-1071-SRG0389

Ø3.5 x 13 mm

70-1071-SRG0390

Ø4.3/5.0 x 8 mm


70-1071-SRG0391

Ø4.3/5.0 x 10 mm

70-1071-SRG0392

Ø4.3/5.0 x 13 mm

70-1071-SRG0393

	Description	Code
Glidewell HT™ Implant Guided Shaping Drill		
 70-1071-SRG0326	Ø3.0 x 11.5 mm	70-1071-SRG0326
 70-1071-SRG0327	Ø3.0 x 13 mm	70-1071-SRG0327
 70-1071-SRG0328	Ø3.0 x 16 mm	70-1071-SRG0328
 70-1071-SRG0329	Ø3.5 x 8 mm	70-1071-SRG0329
 70-1071-SRG0330	Ø3.5 x 10 mm	70-1071-SRG0330
 70-1071-SRG0331	Ø3.5 x 11.5 mm	70-1071-SRG0331
 70-1071-SRG0332	Ø3.5 x 13 mm	70-1071-SRG0332
 70-1071-SRG0333	Ø3.5 x 16 mm	70-1071-SRG0333
 70-1071-SRG0334	Ø4.3 x 8 mm	70-1071-SRG0334
 70-1071-SRG0335	Ø4.3 x 10 mm	70-1071-SRG0335
 70-1071-SRG0336	Ø4.3 x 11.5 mm	70-1071-SRG0336
 70-1071-SRG0337	Ø4.3 x 13 mm	70-1071-SRG0337
 70-1071-SRG0338	Ø4.3 x 16 mm	70-1071-SRG0338

Guided Drills



70-1071-SRG0339



70-1071-SRG0340



70-1071-SRG0341



70-1071-SRG0342



70-1071-SRG0343

Description

Code

Glidewell HT™ Implant Guided Shaping Drills

Ø5.0 x 8 mm

70-1071-SRG0339

Ø5.0 x 10 mm

70-1071-SRG0340

Ø5.0 x 11.5 mm

70-1071-SRG0341

Ø5.0 x 13 mm

70-1071-SRG0342

Ø5.0 x 16 mm

70-1071-SRG0343

Guided Taps & Accessories



70-1071-SRG0362



70-1071-SRG0363



70-1071-SRG0364



70-1071-SRG0365

Glidewell HT™ Implant Guided Screw Tap

Ø3.0 mm

70-1071-SRG0362

Ø3.5 mm


70-1071-SRG0363

Ø4.3 mm

70-1071-SRG0364

Ø5.0 mm

70-1071-SRG0365

	Description	Code
Glidewell HT™ Implant Guided Mount		
 70-1071-SRG0366	Ø3.0 Implant	70-1071-SRG0366
 70-1071-SRG0367	Ø3.5 Implant	70-1071-SRG0367
 70-1071-SRG0368	Ø4.3 Implant	70-1071-SRG0368
 70-1071-SRG0369	Ø5.0 Implant	70-1071-SRG0369
<i>(Guided Mounts include captured screw.)</i>		
 70-1071-SRG0307	Glidewell HT™ Implant Surgical Adaptor - Short	70-1071-SRG0307
 70-1071-SRG0306	Glidewell HT™ Implant Surgical Adaptor	70-1071-SRG0306
 Torque Wrench Driver Handle Hex Driver Handle Square 70-C35120	Glidewell HT™ Torque Wrench with Driver Handle, Hex and Square	70-C35120
 70-1071-SRG0374	Glidewell HT™ Implant Mount Wrench	70-1071-SRG0374
 70-1190-PRC0055	Glidewell HT™ Prosthetic Driver - 24 mm	70-1190-PRC0055
 70-1071-SRG0399	Glidewell HT™ Implant Key - Ø4.3/5.0 to Ø3.5	70-1071-SRG0399

Prosthetic Instrumentation Kit

Features and Benefits

The prosthetic kit is designed specifically for restorative doctors, organizing all of the necessary prosthetic tooling to help simplify restorative procedures and instrument sterilization.

All instrumentation is manufactured in the U.S.A., Switzerland or Denmark. For specific country of origin, please refer to the individual product label.

Sterility

Please refer to the Glidewell HT Implant System Restorative Manual for cleaning and sterilization instructions.

Description	Code
Glidewell™ Prosthetic Kit	70-1071-PRS0004



Prosthetic Instrumentation Kit Overview

Glidewell HT™
IMPLANT SYSTEM

11 
Multi-Unit Driver
70-1071-SRG0234


Glidewell IH™
IMPLANT SYSTEM

1 
Prosthetic Driver - 18 mm
70-1190-PRC0054

2 
Prosthetic Driver - 24 mm
70-1190-PRC0055

3 
Prosthetic Driver - 30 mm
70-1190-PRC0056

4 
Prosthetic Driver - 38 mm
70-1190-PRC0057




10 
Prosthetic Driver, .050 Hex - 18 mm
70-1195-PRC0001

9 
Prosthetic Driver, .050 Hex - 24 mm
70-1195-PRC0002

8 
Prosthetic Driver, .050 Hex - 30 mm
70-1195-PRC0003

7 
Prosthetic Driver, .050 Hex - 38 mm
70-1195-PRC0004

5 
Driver Handle, Square
70-C35119

6 
Torque Wrench
70-C35119

Description

Code



Torque Wrench




Driver Handle Square
70-C35119



70-1071-SRG0234



70-1190-PRC0054



70-1190-PRC0055



70-1190-PRC0056



70-1190-PRC0057



70-1195-PRC0001



70-1195-PRC0002



70-1195-PRC0003



70-1195-PRC0004

Glidewell HT™ Torque Wrench with Driver Handle, Square

70-C35119

Newport Surgical™ Multi-Unit Driver, 4 mm Square

70-1071-SRG0234

Glidewell HT™ Prosthetic Driver - 18 mm

70-1190-PRC0054

Glidewell HT™ Prosthetic Driver - 24 mm

70-1190-PRC0055

Glidewell HT™ Prosthetic Driver - 30 mm

70-1190-PRC0056

Glidewell HT™ Prosthetic Driver - 38 mm

70-1190-PRC0057

Glidewell IH™ Prosthetic Driver - 18 mm

70-1195-PRC0001

Glidewell IH™ Prosthetic Driver - 24 mm

70-1195-PRC0002

Glidewell IH™ Prosthetic Driver - 30 mm

70-1195-PRC0003

Glidewell IH™ Prosthetic Driver - 38 mm



70-1195-PRC0004

Product Packaging

Glidewell HT™ Implants are shipped sterile. They should not be resterilized. They are for single use only, prior to the expiration date. Do not use implants if the packaging has been compromised or previously opened. Do not handle implant surfaces directly. Users are advised to visually inspect packaging to ensure seals and contents are intact prior to use. Please refer to the individual product label for all relevant product information and cautions.

Explanation of Label Codes:

1. Official product description
2. Date of Manufacture (YYYY-MM-DD)
3. Catalog Number
4. Lot/Batch Number
5. By Prescription Only
6. Do Not Use if Package is Damaged
7. Sterile with Gamma Radiation
8. Use-by Date
9. Consult Instructions for Use
10. Do not Re-use
11. Do not Resterilize
12. Manufacturer
13. Country of origin
14. Unique Device Identification (UDI)
15. Non-Sterile
16. QR code for IFU website



Ordering Information

In the U.S.: Order at glidewelldirect.com or call Glidewell Direct at 888-303-3975. Our product specialists are committed to answering questions in a timely fashion to ensure your ordering is easy and efficient. We are available Monday—Friday from 6 a.m.—5 p.m. (*Pacific*).

Shipping Policy

- Orders placed after 3 p.m. (*Pacific*) will be processed on the following business day. Business days do not include Saturdays, Sundays, or U.S. holidays.
- Online shopping cart available to U.S. customers only.

Terms

All accounts are payable within 30 days of invoice date. Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance. We accept American Express, Visa, MasterCard, and Discover. All prices are subject to change without notice.

Product Return Policy

Products may be returned at the customer's expense for credit within 30 days of invoice date. All returned products must meet the following conditions:

- A copy of the original invoice must accompany the products.
- Products must be packaged to arrive at the seller's facility undamaged.
- Discontinued, obsolete, expired, damaged, or opened items will not be accepted for return.
- Amount credited will be based on invoice price, less 15 percent for restocking fee.
- Shipping charges are the responsibility of the customer and will not be credited.

Product & Pricing Changes

Because products and equipment are continually undergoing refinement in design and manufacturing methods, we reserve the right to improve, modify, or discontinue products and equipment or change pricing at any time without incurring any obligation and without prior notice.

Warranty

Limited Warranty—Prismatik Dentalcraft, Inc.

Prismatik Dentalcraft, Inc. ("Prismatik"), is the manufacturer of dental products (the "product"), including Glidewell HT™ Implants ("implants"). Prismatik warrants the Glidewell HT Implant for the life of the patient originally receiving the implant from the date of placement, and for a period of six (6) months for ceramic blanks and any other product ("the warranty period"). Prismatik will at its option replace or refund the purchase price of any product, to the original purchaser ("user"), that is returned due to defects in material and manufacture.

NO GUARANTEE OR WARRANTY IS IMPLIED OTHER THAN EXPRESSLY STATED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

Prismatik shall not be liable for any incidental or consequential damages, whether foreseeable or not, caused by defects in the product or dental devices produced using said product. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, user's exclusive remedy and Prismatik's sole obligation shall be replacement or refund of the purchase price of the product. For replacement or refund under this warranty, the original purchaser shall send the product at its own expense, postage prepaid, to the seller.

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
IMPLANTS					
70-1189-IMP0001	Implant Ø3.0 x 11.5 mm	6	70-1071-SRG0415	Shaping Drill Ø6.0 x 8 mm	21
70-1189-IMP0002	Implant Ø3.0 x 13 mm	6	70-1071-SRG0416	Shaping Drill Ø6.0 x 10 mm	21
70-1189-IMP0003	Implant Ø3.0 x 16 mm	6	70-1071-SRG0417	Shaping Drill Ø6.0 x 11.5 mm	21
70-1189-IMP0004	Implant Ø3.5 x 8 mm	6	70-1071-SRG0288	Shaping Drill Ø7.0 x 8 mm	21
70-1189-IMP0005	Implant Ø3.5 x 10 mm	6	70-1071-SRG0289	Shaping Drill Ø7.0 x 10 mm	21
70-1189-IMP0006	Implant Ø3.5 x 11.5 mm	6	70-1071-SRG0290	Shaping Drill Ø7.0 x 11.5 mm	21
70-1189-IMP0007	Implant Ø3.5 x 13 mm	6	70-1071-SRG0264	Twist Drill Ø2.4/1.5 mm	20
70-1189-IMP0008	Implant Ø3.5 x 16 mm	6	70-1071-SRG0266	Twist Drill Ø1.5 mm x 8 mm	20
70-1189-IMP0009	Implant Ø4.3 x 8 mm	6	70-1071-SRG0268	Twist Drill Ø2.8/2.4 mm x 8 mm	20
70-1189-IMP0010	Implant Ø4.3 x 10 mm	6	70-1071-SRG0269	Twist Drill Ø2.8/2.4 mm x 10 mm	20
70-1189-IMP0011	Implant Ø4.3 x 11.5 mm	6	70-1071-SRG0270	Twist Drill Ø2.8/2.4 mm x 11.5 mm	20
70-1189-IMP0012	Implant Ø4.3 x 13 mm	6	70-1071-SRG0271	Twist Drill Ø2.8/2.4 mm x 13 mm	20
70-1189-IMP0013	Implant Ø4.3 x 16 mm	6	70-1071-SRG0272	Twist Drill Ø2.8/2.4 mm x 16 mm	20
70-1189-IMP0014	Implant Ø5.0 x 8 mm	6	INSTRUMENTATION - DRIVERS & ACCESSORIES		
70-1189-IMP0015	Implant Ø5.0 x 10 mm	6	70-C35119	Torque Wrench with Driver Handle, Square	31
70-1189-IMP0016	Implant Ø5.0 x 11.5 mm	6	70-C35120	Torque Wrench with Driver Handle, Hex and Square ...	22, 29
70-1189-IMP0017	Implant Ø5.0 x 13 mm	6	70-1071-SRG0014	Drill Extender	22
70-1189-IMP0018	Implant Ø5.0 x 16 mm	6	70-1071-SRG0296	Driver Ø3.0, Long	22
70-1189-IMP0022	Implant Ø6.0 x 8 mm	6	70-1071-SRG0297	Driver Ø3.5/4.3, Short	22
70-1189-IMP0023	Implant Ø6.0 x 10 mm	6	70-1071-SRG0298	Driver Ø3.5/4.3, Long	22
70-1189-IMP0024	Implant Ø6.0 x 11.5 mm	6	70-1071-SRG0299	Driver Ø5.0, Short	22
70-1189-IMP0019	Implant Ø7.0 x 8 mm	6	70-1071-SRG0300	Driver Ø5.0, Long	22
70-1189-IMP0020	Implant Ø7.0 x 10 mm	6	70-1071-SRG0419	Driver Ø6.0, Short	22
70-1189-IMP0021	Implant Ø7.0 x 11.5 mm	6	70-1071-SRG0301	Driver Ø7.0, Short	21
INSTRUMENTATION			70-1071-SRG0374	Mount Wrench	29
70-1071-SRG0395	Guided Surgical Kit	23	70-1071-SRG0308	Direction Indicator Ø3.0 mm	22
70-1071-PRS0004	Glidewell Prosthetic Kit	30	70-1071-SRG0309	Direction Indicator Ø3.5 mm	22
70-1071-SRG0396	Surgical Kit	18	70-1071-SRG0310	Direction Indicator Ø4.3 mm	22
INSTRUMENTATION - DRILLS & TAPS			70-1071-SRG0311	Direction Indicator Ø5.0 mm	22
70-1071-SRG0291	Screw Tap Ø3.0 mm	21	70-1071-SRG0413	Direction Indicator Ø6.0 mm	22
70-1071-SRG0292	Screw Tap Ø3.5 mm	21	70-1071-SRG0312	Direction Indicator Ø7.0 mm	22
70-1071-SRG0293	Screw Tap Ø4.3 mm	21	70-1190-PRC0054	Glidewell HT™ Prosthetic Driver - 18 mm	22, 31
70-1071-SRG0294	Screw Tap Ø5.0 mm	21	70-1190-PRC0055	Glidewell HT™ Prosthetic Driver - 24 mm	29, 31
70-1071-SRG0414	Screw Tap Ø6.0 mm	21	70-1190-PRC0056	Glidewell HT™ Prosthetic Driver - 30 mm	22, 31
70-1071-SRG0295	Screw Tap Ø7.0 mm	21	70-1190-PRC0057	Glidewell HT™ Prosthetic Driver - 38 mm	31
70-1071-SRG0273	Shaping Drill Ø3.5 x 8 mm	20	70-1195-PRC0001	Glidewell IH™ Prosthetic Driver - 18 mm	31
70-1071-SRG0274	Shaping Drill Ø3.5 x 10 mm	20	70-1195-PRC0002	Glidewell IH™ Prosthetic Driver - 24 mm	31
70-1071-SRG0275	Shaping Drill Ø3.5 x 11.5 mm	20	70-1195-PRC0003	Glidewell IH™ Prosthetic Driver - 30 mm	31
70-1071-SRG0276	Shaping Drill Ø3.5 x 13 mm	20	70-1195-PRC0004	Glidewell IH™ Prosthetic Driver - 38 mm	31
70-1071-SRG0277	Shaping Drill Ø3.5 x 16 mm	20	70-1071-SRG0234	Newport Surgical Multi-Unit Driver, 4 mm Square	31
70-1071-SRG0278	Shaping Drill Ø4.3 x 8 mm	20	70-1071-RAD0001	Radiographic Template	22
70-1071-SRG0279	Shaping Drill Ø4.3 x 10 mm	20	70-1071-SRG0306	Surgical Adaptor	22, 29
70-1071-SRG0280	Shaping Drill Ø4.3 x 11.5 mm	20	70-1071-SRG0307	Surgical Adaptor, Short	29
70-1071-SRG0281	Shaping Drill Ø4.3 x 13 mm	20	70-1071-SRG0399	Implant Key Ø4.3/5.0 to Ø3.5	29
70-1071-SRG0282	Shaping Drill Ø4.3 x 16 mm	20	INSTRUMENTATION - GUIDED SURGERY		
70-1071-SRG0283	Shaping Drill Ø5.0 x 8 mm	20	70-1071-SRG0317	Guided Alignment Drill - Ø3.0 mm	26
70-1071-SRG0284	Shaping Drill Ø5.0 x 10 mm	20	70-1071-SRG0318	Guided Alignment Drill - Ø3.5 mm	26
70-1071-SRG0285	Shaping Drill Ø5.0 x 11.5 mm	21	70-1071-SRG0319	Guided Alignment Drill - Ø4.3/5.0 mm	26
70-1071-SRG0286	Shaping Drill Ø5.0 x 13 mm	21	70-1071-SRG0366	Guided Mount - Ø3.0 Implant	29
70-1071-SRG0287	Shaping Drill Ø5.0 x 16 mm	21	70-1071-SRG0367	Guided Mount - Ø3.5 Implant	29
			70-1071-SRG0368	Guided Mount - Ø4.3 Implant	29

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
70-1071-SRG0369.....	Guided Mount - Ø5.0 Implant.....	29	70-1190-PRA0021.....	Non-Engaging Temporary Abutment - Ø3.0 Implant.....	8
70-1071-SRG0387.....	Guided Pilot Drill - Ø3.0 x 8 mm.....	26	70-1190-PRA0022.....	Non-Engaging Temporary Abutment - Ø3.5/4.3 Implant.....	8
70-1071-SRG0388.....	Guided Pilot Drill - Ø3.5 x 8 mm.....	26	70-1190-PRA0023.....	Non-Engaging Temporary Abutment - Ø5.0 Implant.....	8
70-1071-SRG0389.....	Guided Pilot Drill - Ø3.5 x 10 mm.....	26	70-1190-PRA0063.....	Non-Engaging Temporary Abutment - Ø6.0 Implant.....	8
70-1071-SRG0390.....	Guided Pilot Drill - Ø3.5 x 13 mm.....	26	70-1190-PRA0024.....	Non-Engaging Temporary Abutment - Ø7.0 Implant.....	8
70-1071-SRG0391.....	Guided Pilot Drill - Ø4.3/5.0 x 8 mm.....	26	70-1190-PRA0025.....	Temporary Abutment - Ø3.0 Implant.....	7
70-1071-SRG0392.....	Guided Pilot Drill - Ø4.3/5.0 x 10 mm.....	26	70-1190-PRA0026.....	Temporary Abutment - Ø3.5/4.3 Implant.....	7
70-1071-SRG0393.....	Guided Pilot Drill - Ø4.3/5.0 x 13 mm.....	26	70-1190-PRA0027.....	Temporary Abutment - Ø5.0 Implant.....	7
70-1071-SRG0362.....	Guided Screw Tap - Ø3.0 mm.....	28	70-1190-PRA0062.....	Temporary Abutment - Ø6.0 Implant.....	7
70-1071-SRG0363.....	Guided Screw Tap - Ø3.5 mm.....	28	70-1190-PRA0028.....	Temporary Abutment - Ø7.0 Implant.....	7
70-1071-SRG0364.....	Guided Screw Tap - Ø4.3 mm.....	28			
70-1071-SRG0365.....	Guided Screw Tap - Ø5.0 mm.....	27			
70-1071-SRG0326.....	Guided Shaping Drill - Ø3.0 x 11.5 mm.....	27			
70-1071-SRG0327.....	Guided Shaping Drill - Ø3.0 x 13 mm.....	27			
70-1071-SRG0328.....	Guided Shaping Drill - Ø3.0 x 16 mm.....	27			
70-1071-SRG0329.....	Guided Shaping Drill - Ø3.5 x 8 mm.....	27			
70-1071-SRG0330.....	Guided Shaping Drill - Ø3.5 x 10 mm.....	27			
70-1071-SRG0331.....	Guided Shaping Drill - Ø3.5 x 11.5 mm.....	27			
70-1071-SRG0332.....	Guided Shaping Drill - Ø3.5 x 13 mm.....	27			
70-1071-SRG0333.....	Guided Shaping Drill - Ø3.5 x 16 mm.....	27			
70-1071-SRG0334.....	Guided Shaping Drill - Ø4.3 x 8 mm.....	27			
70-1071-SRG0335.....	Guided Shaping Drill - Ø4.3 x 10 mm.....	27			
70-1071-SRG0336.....	Guided Shaping Drill - Ø4.3 x 11.5 mm.....	27			
70-1071-SRG0337.....	Guided Shaping Drill - Ø4.3 x 13 mm.....	27			
70-1071-SRG0338.....	Guided Shaping Drill - Ø4.3 x 16 mm.....	27			
70-1071-SRG0339.....	Guided Shaping Drill - Ø5.0 x 8 mm.....	28			
70-1071-SRG0340.....	Guided Shaping Drill - Ø5.0 x 10 mm.....	28			
70-1071-SRG0341.....	Guided Shaping Drill - Ø5.0 x 11.5 mm.....	28			
70-1071-SRG0342.....	Guided Shaping Drill - Ø5.0 x 13 mm.....	28			
70-1071-SRG0343.....	Guided Shaping Drill - Ø5.0 x 16 mm.....	28			
70-1071-SRG0313.....	Guided Tissue Punch - Ø3.0 mm.....	26			
70-1071-SRG0314.....	Guided Tissue Punch - Ø3.5 mm.....	26			
70-1071-SRG0315.....	Guided Tissue Punch - Ø4.3 mm.....	26			
70-1071-SRG0316.....	Guided Tissue Punch - Ø5.0 mm.....	26			
PROSTHETICS - HEALING/TEMPORIZATION			PROSTHETICS - IMPRESSIONS/ANALOGS		
70-1190-IMP0001.....	Healing Abutment 3 mmH - Ø3.0 Implant.....	7	70-1190-PRC0035.....	Analog - Ø3.0 Implant.....	10
70-1190-IMP0002.....	Healing Abutment 5 mmH - Ø3.0 Implant.....	7	70-1190-PRC0036.....	Analog - Ø3.5 Implant.....	10
70-1190-IMP0003.....	Healing Abutment 7 mmH - Ø3.0 Implant.....	7	70-1190-PRC0037.....	Analog - Ø4.3 Implant.....	10
70-1190-IMP0004.....	Healing Abutment 3 mmH - Ø3.5 Implant.....	7	70-1190-PRC0038.....	Analog - Ø5.0 Implant.....	10
70-1190-IMP0005.....	Healing Abutment 5 mmH - Ø3.5 Implant.....	7	70-1190-PRC0062.....	Analog - Ø6.0 Implant.....	10
70-1190-IMP0006.....	Healing Abutment 7 mmH - Ø3.5 Implant.....	7	70-1190-PRC0039.....	Analog - Ø7.0 Implant.....	10
70-1190-IMP0007.....	Healing Abutment 3 mmH - Ø4.3 Implant.....	7	70-1190-PRC0001.....	Closed-Tray Impression Coping 3 mmH - Ø3.0.....	8
70-1190-IMP0008.....	Healing Abutment 5 mmH - Ø4.3 Implant.....	7	70-1190-PRC0002.....	Closed-Tray Impression Coping 5 mmH - Ø3.0.....	8
70-1190-IMP0009.....	Healing Abutment 7 mmH - Ø4.3 Implant.....	7	70-1190-PRC0003.....	Closed-Tray Impression Coping 7 mmH - Ø3.0.....	8
70-1190-IMP0010.....	Healing Abutment 3 mmH - Ø5.0 Implant.....	7	70-1190-PRC0004.....	Closed-Tray Impression Coping 3 mmH - Ø3.5.....	8
70-1190-IMP0011.....	Healing Abutment 5 mmH - Ø5.0 Implant.....	7	70-1190-PRC0005.....	Closed-Tray Impression Coping 5 mmH - Ø3.5.....	8
70-1190-IMP0012.....	Healing Abutment 7 mmH - Ø5.0 Implant.....	7	70-1190-PRC0006.....	Closed-Tray Impression Coping 7 mmH - Ø3.5.....	8
70-1190-IMP0015.....	Healing Abutment 3 mmH - Ø6.0 Implant.....	7	70-1190-PRC0007.....	Closed-Tray Impression Coping 3 mmH - Ø4.3.....	8
70-1190-IMP0016.....	Healing Abutment 5 mmH - Ø6.0 Implant.....	7	70-1190-PRC0008.....	Closed-Tray Impression Coping 5 mmH - Ø4.3.....	8
70-1190-IMP0013.....	Healing Abutment 3 mmH - Ø7.0 Implant.....	7	70-1190-PRC0009.....	Closed-Tray Impression Coping 7 mmH - Ø4.3.....	8
70-1190-IMP0014.....	Healing Abutment 5 mmH - Ø7.0 Implant.....	7	70-1190-PRC0010.....	Closed-Tray Impression Coping 3 mmH - Ø5.0.....	8
			70-1190-PRC0011.....	Closed-Tray Impression Coping 5 mmH - Ø5.0.....	8
			70-1190-PRC0012.....	Closed-Tray Impression Coping 7 mmH - Ø5.0.....	8
			70-1190-PRC0058.....	Closed-Tray Impression Coping 3 mmH - Ø6.0.....	8
			70-1190-PRC0059.....	Closed-Tray Impression Coping 5 mmH - Ø6.0.....	8
			70-1190-PRC0013.....	Closed-Tray Impression Coping 3 mmH - Ø7.0.....	8
			70-1190-PRC0014.....	Closed-Tray Impression Coping 5 mmH - Ø7.0.....	8
			70-1190-PRC0015.....	Open-Tray Impression Coping 3 mmH - Ø3.0.....	9
			70-1190-PRC0016.....	Open-Tray Impression Coping 5 mmH - Ø3.0.....	9
			70-1190-PRC0017.....	Open-Tray Impression Coping 7 mmH - Ø3.0.....	9
			70-1190-PRC0018.....	Open-Tray Impression Coping 3 mmH - Ø3.5.....	9
			70-1190-PRC0019.....	Open-Tray Impression Coping 5 mmH - Ø3.5.....	9
			70-1190-PRC0020.....	Open-Tray Impression Coping 7 mmH - Ø3.5.....	9
			70-1190-PRC0021.....	Open-Tray Impression Coping 3 mmH - Ø4.3.....	9
			70-1190-PRC0022.....	Open-Tray Impression Coping 5 mmH - Ø4.3.....	9
			70-1190-PRC0023.....	Open-Tray Impression Coping 7 mmH - Ø4.3.....	9
			70-1190-PRC0024.....	Open-Tray Impression Coping 3 mmH - Ø5.0.....	9
			70-1190-PRC0025.....	Open-Tray Impression Coping 5 mmH - Ø5.0.....	9
			70-1190-PRC0026.....	Open-Tray Impression Coping 7 mmH - Ø5.0.....	9
			70-1190-PRC0060.....	Open-Tray Impression Coping 3 mmH - Ø6.0.....	9
			70-1190-PRC0061.....	Open-Tray Impression Coping 5 mmH - Ø6.0.....	9
			70-1190-PRC0027.....	Open-Tray Impression Coping 3 mmH - Ø7.0.....	9
			70-1190-PRC0028.....	Open-Tray Impression Coping 5 mmH - Ø7.0.....	9
			70-1190-PRA0054.....	Titanium Scan Body - Ø3.0 Implant.....	10
			70-1190-PRA0055.....	Titanium Scan Body - Ø3.5/4.3 Implant.....	10

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
70-1190-PRA0056	Titanium Scan Body - Ø5.0 Implant.....	10	70-1190-PRS0020	Multi-Unit Abutment 2.5 mmH - Ø6.0 Implant.....	13
70-1190-PRA0064	Titanium Scan Body - Ø6.0 Implant.....	10	70-1190-PRS0021	Multi-Unit Abutment 3.5 mmH - Ø6.0 Implant.....	13
70-1190-PRA0057	Titanium Scan Body - Ø7.0 Implant.....	10	70-1190-PRS0022	Multi-Unit Abutment 4.5 mmH - Ø6.0 Implant.....	13
PROSTHETICS - LOCATOR ABUTMENT SYSTEM			70-1190-PRS0018	Multi-Unit Coping 4.5 mmH, Non-Engaging.....	14
70-1075-ZST0001	35 Ncm Torque Wrench Insert Driver (15 mm).....	17	70-1074-PRS0016	Multi-Unit Abutment Analog.....	15
70-1075-ZST0002	35 Ncm Torque Wrench Insert Driver (21 mm).....	17	70-1071-SRG0234	Newport Surgical Multi-Unit Driver, 4mm Square.....	15
70-1075-ZST0026	Angle Measurement Guide.....	17	70-1074-PRS0018	Multi-Unit Gold Bar Coping.....	15
70-1075-ZST0004	Core Tool (Complete).....	17	70-1074-PRS0030	Multi-Unit Gold/Plastic Coping.....	15
70-1075-ZST0003	Core Tool Abutment Driver (Gold End).....	17	70-1074-PRS0019	Multi-Unit Guide Pin.....	15
70-1075-ZST0005	Core Tool Abutment Holder Sleeve (4/pk).....	17	70-1074-PRS0021	Multi-Unit Impression Coping, Closed-Tray.....	15
70-1075-ZST0006	Core Tool Male Removal (Tip End).....	17	70-1074-PRS0022	Multi-Unit Impression Coping, Open-Tray.....	15
70-1075-ZST0025	Extended Range Male Processing Package (2/pk).....	17	70-1074-PRS0069	Multi-Unit Plastic Coping.....	15
70-1075-ZST0009	Female Analog 4 mmD (4/pk).....	17	70-1072-PRC0043	Multi-Unit Prosthetic Screw.....	15
70-1075-ZST0007	Impression Coping (4/pk).....	17	70-1074-PRS0044	Multi-Unit Temporary Healing Cap.....	15
70-1075-ZST0020	Latch Type Torque Wrench Insert Driver (23 mm).....	17	70-1190-PRC0046	Multi-Unit Titanium Scan Body with Prosthetic Screw...	15
70-1075-ZST0021	Latch Type Torque Wrench Insert Driver (29 mm).....	17	70-1074-PRS0057	Multi-Unit Titanium Temporary.....	15
70-1075-ZST0065	Locator Abutment Ø3.0 x 1.0 mmH Cuff.....	16	PROSTHETICS - SCREWS & GUIDE PINS		
70-1075-ZST0074	Locator Abutment Ø3.0 x 2.0 mmH Cuff.....	16	70-1190-PRC0047	Guide Pin - Ø3.0 Implant.....	13
70-1075-ZST0066	Locator Abutment Ø3.0 x 3.0 mmH Cuff.....	16	70-1190-PRC0048	Guide Pin - Ø3.5/4.3 Implant.....	13
70-1075-ZST0032	Locator Abutment Ø3.5/4.3 x 1.0 mmH Cuff.....	16	70-1190-PRC0065	Guide Pin - Ø5.0/6.0 Implant.....	13
70-1075-ZST0034	Locator Abutment Ø3.5/4.3 x 3.0 mmH Cuff.....	16	70-1190-PRC0049	Guide Pin - Ø5.0/7.0 Implant.....	13
70-1075-ZST0036	Locator Abutment Ø3.5/4.3 x 5.0 mmH Cuff.....	16	70-1190-PRC0029	Titanium Screw - Ø3.0 Implant.....	13
70-1075-ZST0071	Locator Abutment Ø5.0 x 1.0 mmH Cuff.....	16	70-1190-PRC0030	Titanium Screw - Ø3.5/4.3 Implant.....	13
70-1075-ZST0038	Locator Abutment Ø5.0 x 3.0 mmH Cuff.....	16	70-1190-PRC0064	Titanium Screw - Ø5.0/6.0 Implant.....	13
70-1075-ZST0040	Locator Abutment Ø5.0 x 5.0 mmH Cuff.....	16	70-1190-PRC0031	Titanium Screw - Ø5.0/7.0 Implant.....	13
70-1075-ZST0024	Male Processing Package (2/pk).....	17	PROSTHETICS - TITANIUM ABUTMENTS		
70-1075-ZST0008	Parallel Post (4/pk).....	17	70-1190-PRA0001	Non-Engaging Ti. Abutment 4.5 mmH - Ø3.0 Implant....	11
70-1075-ZST0010	Processing Spacer (4/pk).....	17	70-1190-PRA0002	Non-Engaging Ti. Abutment 4.5 mmH - Ø3.5/4.3 Implant	11
70-1075-ZST0022	Square Drive Torque Wrench Insert Driver (15 mm).....	17	70-1190-PRA0003	Non-Engaging Ti. Abutment 4.5 mmH - Ø5.0 Implant....	11
70-1075-ZST0023	Square Drive Torque Wrench Insert Driver (21 mm).....	17	70-1190-PRA0065	Non-Engaging Ti. Abutment 4.5 mmH - Ø6.0 Implant....	11
PROSTHETICS - MULTI-UNIT ABUTMENT SYSTEM			70-1190-PRA0004	Non-Engaging Ti. Abutment 4.5 mmH - Ø7.0 Implant.....	11
70-1190-PRS0001	17° Multi-Unit Abutment 2.5 mmH - Ø3.5/4.3 Implant....	14	70-1190-PRA0005	Titanium Abutment 4.5 mmH - Ø3.0 Implant.....	11
70-1190-PRS0002	17° Multi-Unit Abutment 2.5 mmH - Ø5.0 Implant.....	14	70-1190-PRA0006	Titanium Abutment 4.5 mmH - Ø3.5/4.3 Implant.....	11
70-1190-PRS0003	17° Multi-Unit Abutment 3.5 mmH - Ø3.5/4.3 Implant....	14	70-1190-PRA0007	Titanium Abutment 4.5 mmH - Ø5.0 Implant.....	11
70-1190-PRS0004	17° Multi-Unit Abutment 3.5 mmH - Ø5.0 Implant.....	14	70-1190-PRA0066	Titanium Abutment 4.5 mmH - Ø6.0 Implant.....	11
70-1190-PRS0005	30° Multi-Unit Abutment 3.5 mmH - Ø3.5/4.3 Implant....	14	70-1190-PRA0008	Titanium Abutment 4.5 mmH - Ø7.0 Implant.....	11
70-1190-PRS0006	30° Multi-Unit Abutment 3.5 mmH - Ø5.0 Implant.....	14	70-1190-PRA0009	Titanium Abutment 6 mmH - Ø3.0 Implant.....	11
70-1190-PRS0007	30° Multi-Unit Abutment 4.5 mmH - Ø3.5/4.3 Implant....	14	70-1190-PRA0010	Titanium Abutment 6 mmH - Ø3.5/4.3 Implant.....	11
70-1190-PRS0008	30° Multi-Unit Abutment 4.5 mmH - Ø5.0 Implant.....	14	70-1190-PRA0011	Titanium Abutment 6 mmH - Ø5.0 Implant.....	11
70-1190-PRC0032	Angled Multi-Unit Abutment Screw - Ø3.5/4.3 Implant..	14	70-1190-PRA0067	Titanium Abutment 6 mmH - Ø6.0 Implant.....	11
70-1190-PRC0033	Angled Multi-Unit Abutment Screw - Ø5.0 Implant.....	14	70-1190-PRA0012	Titanium Abutment 6 mmH - Ø7.0 Implant.....	11
70-1190-PRS0009	Multi-Unit Abutment 1.5 mmH - Ø3.5/4.3 Implant.....	13	PROSTHETICS - TITANIUM ESTHETIC ABUTMENTS		
70-1190-PRS0010	Multi-Unit Abutment 1.5 mmH - Ø5.0 Implant.....	13	70-1190-PRA0042	15° Anterior Esthetic Abutment - Ø3.0 Implant.....	11
70-1190-PRS0011	Multi-Unit Abutment 2.5 mmH - Ø3.0 Implant.....	13	70-1190-PRA0043	15° Anterior Esthetic Abutment - Ø3.5/4.3 Implant.....	11
70-1190-PRS0012	Multi-Unit Abutment 2.5 mmH - Ø3.5/4.3 Implant.....	13	70-1190-PRA0044	15° Anterior Esthetic Abutment - Ø5.0 Implant.....	11
70-1190-PRS0013	Multi-Unit Abutment 2.5 mmH - Ø5.0 Implant.....	13	70-1190-PRA0045	15° Posterior Esthetic Abutment - Ø5.0 Implant.....	11
70-1190-PRS0014	Multi-Unit Abutment 3.5 mmH - Ø3.0 Implant.....	13	70-1190-PRA0037	Anterior Esthetic Abutment - Ø3.0 Implant.....	10
70-1190-PRS0015	Multi-Unit Abutment 3.5 mmH - Ø3.5/4.3 Implant.....	13	70-1190-PRA0038	Anterior Esthetic Abutment - Ø3.5/4.3 Implant.....	10
70-1190-PRS0016	Multi-Unit Abutment 3.5 mmH - Ø5.0 Implant.....	13	70-1190-PRA0039	Anterior Esthetic Abutment - Ø5.0 Implant.....	10
70-1190-PRS0017	Multi-Unit Abutment 4.5 mmH - Ø5.0 Implant.....	13	70-1190-PRA0040	Posterior Esthetic Abutment - Ø5.0 Implant.....	10
70-1190-PRS0019	Multi-Unit Abutment 1.5 mmH - Ø6.0 Implant.....	13	70-1190-PRA0075	Posterior Esthetic Abutment - Ø6.0 Implant.....	10

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
70-1190-PRA0041	Posterior Esthetic Abutment - Ø7.0 Implant	10			
PROSTHETICS - UCLA ABUTMENTS					
70-1190-PRA0050	Non-Engaging UCLA Gold - Ø3.0	12			
70-1190-PRA0051	Non-Engaging UCLA Gold - Ø3.5/4.3	12			
70-1190-PRA0052	Non-Engaging UCLA Gold - Ø5.0	12			
70-1190-PRA0073	Non-Engaging UCLA Gold - Ø6.0	12			
70-1190-PRA0053	Non-Engaging UCLA Gold - Ø7.0	12			
70-1190-PRA0033	Non-Engaging UCLA Plastic - Ø3.0	12			
70-1190-PRA0034	Non-Engaging UCLA Plastic - Ø3.5/4.3	12			
70-1190-PRA0035	Non-Engaging UCLA Plastic - Ø5.0	12			
70-1190-PRA0071	Non-Engaging UCLA Plastic - Ø6.0	12			
70-1190-PRA0036	Non-Engaging UCLA Plastic - Ø7.0	12			
70-1190-PRA0046	UCLA Gold - Ø3.0	12			
70-1190-PRA0047	UCLA Gold - Ø3.5/4.3	12			
70-1190-PRA0048	UCLA Gold - Ø5.0	12			
70-1190-PRA0070	UCLA Gold - Ø6.0	12			
70-1190-PRA0049	UCLA Gold - Ø7.0	12			
70-1190-PRA0029	UCLA Plastic - Ø3.0	12			
70-1190-PRA0030	UCLA Plastic - Ø3.5/4.3	12			
70-1190-PRA0031	UCLA Plastic - Ø5.0	12			
70-1190-PRA0074	UCLA Plastic - Ø6.0	12			
70-1190-PRA0032	UCLA Plastic - Ø7.0	12			

1. Kerr M, Allen B, Park N. Clinical and radiographic evaluation of tapered implants with an aggressive reverse buttress thread and crestal microthreads: a retrospective study. For the full report, visit glidewell.com/ht-2-year.
2. Gracis S, Llobell A, Bichacho N, Jahangiri L, Ferencz JL. The Influence of Implant Neck Features and Abutment Diameter on Hard and Soft Tissues Around Single Implants Placed in Healed Ridges: Clinical Criteria for Selection. *Int J Periodontics Restorative Dent*. 2020 Jan/Feb;40(1):39-48.
3. Germanier Y, Tosatti S, Brogini N, Textor M, Buser D. Enhanced bone apposition around biofunctionalized sandblasted and acid-etched titanium implant surfaces. *Clin Oral Implants Res*. 2006 Jun;17(3):251-7.
4. Dohan Ehrenfest DM, Coelho PG, Kang BS, Sul YT, Albrektsson T. Classification of osseointegrated implant surfaces: materials, chemistry and topography. *Trends Biotechnol*. 2010 Apr;28(4):198-206. doi: 10.1016/j.tibtech.2009.12.003. Epub 2010 Jan 29.
5. Kim YK, Kim SG, Kim JH, Yi YJ, Yun PY. Prospective study of tapered resorbable blasting media surface implant stability in the maxillary posterior area. *Oral Surg Oral Med Oral Pathol Oral Radiol*. 2012 Jul;114(1):e19-24. doi: 10.1016/j.tripleo.2011.08.028. Epub 2012 Feb 28.
6. Abuhussein H, Pagni G, Rebaudi A, Wang HL. The effect of thread pattern upon implant osseointegration. *Clin Oral Implants Res*. 2010 Feb;21(2):129-36. doi: 10.1111/j.1600-0501.2009.01800.x. Epub 2009 Aug 25.
7. Hermann JS, Jones AA, Bakaeen LG, Buser D, Schoolfield JD, Cochran DL. Influence of a machined collar on crestal bone changes around titanium implants: a histometric study in the canine mandible. *J Periodontol*. 2011;82:1329-1338.
8. Ormianer Z, Palti A. Retrospective clinical evaluation of tapered screw-vent implants: results after up to eight years of clinical function. *J Oral Implantol*. 2008;34(3):150-60.
9. Goswami M. Comparison of crestal bone loss among two implant crest module designs. *MJAFI*. 2009;65:319-322.
10. Torroella-Saura G, Mareque-Bueno J, Cabratosa-Termes J, Hernández-Alfaro F, Ferrés-Padró E, Calvo-Guirado JL. Effect of implant design in immediate loading. A randomized, controlled, split-mouth, prospective clinical trial. *Clin Oral Implants Res*. 2015 Mar;26(3):240-4.
11. Eraslan O, Inan O. The effect of thread design on stress distribution in a solid screw implant: a 3D finite element analysis. *Clin Oral Investig*. 2010 Aug;14(4):411-6.
12. Alves CC, Neves M. Tapered implants: from indications to advantages. *Int J Periodontics Restorative Dent*. 2009 Apr;29(2):161-7.
13. Moergel M, Rocha S, Messias A, Nicolau P, Guerra F, Wagner W. Radiographic evaluation of conical tapered platform-switched implants in the posterior mandible: 1-year results of a two-center prospective study. *Clin Oral Implants Res*. 2015 Jun 21.

Glidewell HT™

IMPLANT SYSTEM



888-303-3975



glidewell-ht.com



PRISMATIK
DENTALCRAFT, INC.

2144 Michelson Drive • Irvine, CA 92612, USA

